

ULTRASONIC OSTEOTOMY SURGICAL SYSTEM

olea® XD880A




Ultrasonic Osteotome Expert


olea® XD880A is an ultrasonic osteotome that enables safer and more controlled bone removal while sparing soft tissue and neural structures.



SMTP



Drill




olea® XD880A is a valuable tool in cutting, drilling and grinding in spine surgery, neurosurgery, orthopedic surgery, joint surgery and endoscopic surgery.

Non-rotational longitudinal micro-vibration provides for excellent bone removal




Safety selection of Tips

No.1 Bilateral Scalpel with Serrations




No.2 Curette with teeth




No.3 Spiral-shaped Scalpel



No.4 V-shaped Scalpel



No.5 Hook-shaped Scalpel with Serrations




A wide range of advanced tips are available

Safety for Patients


Better protecting for soft tissue

" Low Amplitude + Blunt Tip " efficiently protects the spinal dura mater.




Higher protection

Class I type BF protects against electric shock, which provides a higher electrical safety and protection for doctors and patients.




Appropriate for Doctors

Engineered for
endoscopic, MED and open spine surgery



Efficient cutting , drilling and grinding
with lower noise

A. Customized shape, length and angle of TIP,
fulfill your surgical requirements




B. Designed for inward and outward resection


Inward resection = Outside → inside

Outward resection = Inside → outside


Wide Range of Clinical Applications





Spine surgery





Neurosurgery




Joint surgery





Hand and foot surgery



Pediatric orthopedic surgery




Otolaryngologic surgery



Bone tumor surgery



Maxillofacial and plastic surgery



CONTACT US

Add: Room 2301, Sector B Building 1, No. 29 Qingfeng West Road,
Daxing Bio-medicine Industry Park, Daxing District, Beijing (102629), China
Tel: 86-10-88571203
Fax: 86-10-88572898-8005
Email: marketing@smtpmc.com
Web: www.smtpmc.com

Distributor

Specifications subject to change without
prior notice 2022 SMTP Technology Co ,
Ltd . All rights reserved