


SMT
TECHNOLOGY

agile 

SMT Ultrasonic Surgical System



• wing5™ Series Shears •



Brand New Shears

SAFE HEMOSTASIS GUARANTEE



Unique ACT Technology

EFFICIENT VASCULAR COAGULATION



Original Design

PRODUCED BY NATIONAL INDUSTRY
STANDARD DRAFTER



PRECISE AND DELICATED
BLADE



ERGONOMIC BUTTON
DESIGN

wing5™ Series Shears Scalpels

- Optimize sealing and dissection for vessels, Safe and Reliable
- 24cm/34cm/43cm Shafts available, meeting wide range of surgical requirements
- Faster cutting speed, Lower maximum blade temperature, Minimal thermal injury
- Less smoke providing greater vision in endoscopic surgery
- Less tissue coking ensuring better patient recovery



agile™ Ultrasonic Surgical Generator

- Minimalist design, agile and easy transported, multiple placement options in OR
- Auto intelligent self-testing ensures smooth surgical operation
- ACT™ Adapted Coagulating Technology providing confident vessel sealing
- Equipped with Empty Clip Prompt function, reminding surgeon of safer surgical operations
- Original design, unique 52kHz output frequency, more efficient and faster cutting and coagulation

